

1. Objective

This protocol describes an assay that studies antibody transcytosis across a blood-brain barrier on-a-chip in the MIMETAS OrganoPlate® 3-lane (as [published](#) by Wevers et al., DOI: 10.1186/S12987-018-0108-3).

2. Background

Therapeutic antibodies can be targeted to the brain to treat brain diseases. In this assay, a combination of a target antibody (which is supposed to cross the BBB) and a control antibody (which is not supposed to cross the BBB) is perfused through the apical compartment of the BBB model (see figure 1). After incubation, samples are taken from the basal compartment. The antibody contents of these samples can be analyzed using ELISA or a Mesoscale Discovery system.




Figure 1. The OrganoPlate® 3-lane comprises 40 chips, which are used to culture 40 miniaturized BBB models. The BBB model comprises an endothelial microvessel, grown against an extracellular matrix (ECM) gel, and supporting astrocytes and pericytes. Antibody is perfused through the lumen of the endothelial microvessel (the apical compartment). Samples are then taken from the basal compartment and antibody contents are analyzed.

3. Materials

- OrganoPlate® 3-lane with BBB co-cultures (see figure 2)
 - The culture has an endothelial tubule in the first channel (= apical compartment)
 - The culture has astrocytes and pericytes in the third channel (= basal compartment)
- Medium (i.e. ScienCell™ Astrocyte medium, ScienCell #1801)
- Antibodies
 - Target antibody
 - Control antibody
- Multichannel pipette
- Multichannel tips
- OrganoFlow®




Figure 2. BBB co-culture of endothelial cells (red), astrocytes (green), and pericytes (blue) in a chip of an OrganoPlate® 3-lane.

4. Assay


Prepare antibody solutions

1. Warm medium in a 37°C water bath
2. Prepare the antibody solution in pre-warmed medium at the appropriate concentration
 - a. Prepare 70 µL of antibody solution per chip
 - b. Optimal antibody concentrations depend on the antibody used
 - i. For MEM-189 and anti-HEL antibodies, we used a concentration of 1.25 µM
 - c. Ideally, prepare a solution that contains both a target antibody (i.e. from human) and a control antibody (i.e. from mouse). For example, prepare 1.25 µM target antibody + 1.25 µM control antibody in pre-warmed medium. This way, each target has a control measured from the same chip. If the target and control antibody originate from the same species, you can use different chips or adapt the analysis method.

Wash chips with medium


During the culture of the BBB on-a-chip model, only the top channel is perfused with medium. To assure proper flow profiles and correct readouts, it is essential to first perform a “washing step” that wets all the inlets and outlets.

3. Aspirate media from top inlets and outlets
4. Add 25 µL **medium** to all inlets and outlets (see image on the right)
5. Place plate under an angle and perfuse for 5 minutes
6. Aspirate all media
7. Plate is now ready for starting antibody transcytosis assay



Antibody transcytosis

8. Add 20 µL of **medium** without antibody to the middle inlets and outlets and the bottom inlets and outlets of each chip
9. Add 35 µL of **antibody solution** to the top inlets and outlets
10. Place plate on OrganoFlow® in the incubator for 1 hour
 - a. For different endothelial cell types and different antibodies, the incubation time may need optimization
11. Sample from apical and basal compartments by collecting the contents of the associated wells
 - a. In case you want to sample from one chip at several time points, you can remove the volume you need from the compartment of interest by pipetting (i.e. 10 µL). You can then decide to not adjust the volumes (resulting in a flow through the system) or to replace the volume taken out with fresh medium (and correct for the change in concentration).



MIMETAS product list

Cat. No.	Product Name
MI-AR-CC-01	OrganoReady® Caco-2
9605-400-B	OrganoPlate® 2-lane
4004-400-B	OrganoPlate® 3-lane 40
6405-400-B	OrganoPlate® 3-lane 64
6401-400-B	OrganoPlate® Graft
MI-OFPR-S	OrganoFlow® S
MI-OFPR-L	OrganoFlow® L
MI-OT-1	OrganoTEER®

Contact information

For questions, please contact us through the e-mail addresses stated below

Purchasing: order@mimetas.com

Customer service: info@mimetas.com

Technical support: support@mimetas.com

MIMETAS Europe

J.H. Oortweg 19
2333 CH, Leiden
The Netherlands
Phone: +31 (0)85 888 3161

MIMETAS USA

704 Quince Orchard Road
Suite 260, MD 20878
Gaithersburg, USA
+1 (833) 646-3827

MIMETAS Japan

4F Tekko Building,
1-8-2 Marunouchi, Chiyoda-Ku
Tokyo, 100-0005, Japan
+81 3-6870-7235

This protocol is provided 'as is' and without any warranties, express or implied, including any warranty of merchantability or fitness for a particular purpose or assured results, or that the use of the protocol will not infringe any patent, copyright, trademark, or other proprietary rights. This protocol cannot be used for diagnostic purposes or be resold. The use of this protocol is subject to Mimetas' General Terms and Conditions of Delivery, Purchase and Use.

MIMETAS®, OrganoPlate®, OrganoFlow®, OrganoReady®, and OrganoTEER® are registered trademarks of MIMETAS BV. ScienCell™ is a registered trademark of ScienCell Research Laboratories, Inc.