


MIMETAS
OrganoPlate® family



2-lane 96



3-lane 40



3-lane 64



Graft

Functionalities

	96	40	64	64
Tissue culture chips	96	40	64	64
Apical acces tubules	✓	✓	✓	✓
Basolateral acces tubules	—	✓	✓	✓
Perfusion	✓	✓	✓	✓
3D	✓	✓	✓	✓
Permeability/transport	—	✓	✓	✓
Co-culture	✓	✓	✓	✓
Migration/invasion	—	✓	✓	✓
Angiogenesis	—	✓	✓	✓
Vascularization	—	—	—	✓
OrganoTEER®	—	✓	—	—
Specially designed for automation	—	—	✓	✓